


# Handhabung

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
# How to use

Beware! Do not pull the thumb coach forcefully over the fingers! If it gets stuck, hold under cold water and use detergent

ca. 10x ③ & ④

# Thumbcoach

## Frequent texters, listen up!

**It is remarkable how often and at what speed we use our smart phones in today's communication age. A healthy thumb is therefore all the more important.**

Whatsapp, Facebook, SMS, email, pictures - we can control, answer, comment on or record everything with our smart phone using our thumbs.

**We usually don't notice that we are putting our thumbs under high stress until they start to hurt.**

This unnatural, circular movement which we employ when writing a message is unbalanced and can cause a strain or even tendonitis.

So-called 'smart phone thumb' or 'game pad thumb' is a typical symptom which is occurring ever more frequently and many doctors and physiotherapists are sure that it is becoming a widespread illness.

**THUMBCOACH** was developed to avoid such injuries.

You can find the medical publication "SMS thumb" and much more information at: [www.thumb-coach.com](http://www.thumb-coach.com)

## Function and application of the THUMBCOACH

The shape of the **THUMBCOACH** is shaped to the human hand and holds the four fingers of the hand in a stable position. The thumb thereby remains mobile, tightening and relaxing the rubber band.

This is the exact counter-movement to the strain put on the thumb when operating a smart phone or game pad.

In one training session, the rubber band is stretched 10 to 20 times with the left hand and then with the right hand. The exercise is performed three times in a row. With regular use, your thumb will stay fit.

The **THUMBCOACH** is exclusively for prevention.

If you already have symptoms, contact your doctor or physiotherapist. Incorrect use can also have negative effects.

The **THUMBCOACH** is available in sizes S, M & L